MultiSync® MD Series

21.3" grayscale medical-grade LCD monitors for diagnostic applications



Modern medical technology enables complex examinations. But in the end, it is the radiologist who is responsible for fast and accurate diagnostic findings. To support these findings in a soft copy reading environment, flat-panel displays need to comply with the highest visual requirements. Their image quality has a decisive influence on diagnostic accuracy.

NEC Display Solutions is setting new standards with grayscale displays for filmless diagnostics, with technologies that are available for the first time in this market segment. This introduction of medical-grade displays is backed by years of experience in the development of TFT displays for professional use combined with intensive research done in the medical market.

MultiSync MD Series grayscale displays – in combination with the appropriate graphics controller – are designed for diagnostic imaging and for use in Picture Archiving and Communication System (PACS) applications. Furthermore, diagnostic imaging displays from NEC Di

areas of application and reduced service costs are just a few examples that contribute to cost reduction in healthcare.



Factory calibration and consistency. Using state-of-the-art equipment during production, each NEC MultiSync MD Series monitor is calibrated to the DICOM display function for luminance and to a desired white-

X-Light technology - a view behind the screen

- 1 LCD panel
- 2 Polariser diffusing plate, etc.
- 3 Light guide plate
- 4 Color sensor
- 5 Temperature sensor
- 6 Current data of color sensor and temperature sensor
- 7 Control of 3 inverters
- 8 Calibration of whitepoint via pre-calibration in factory, software (GammaComp MD) or reference value of external sensor
- 9 Light sources, white light using a mixture of different phosphor colors

point (color of the white) -whether for clear base or blue base. NEC Display Solutions' patent-pending X-Light[™] technology (Fig. 1) has the unique ability to control and adjust the luminance and whitepoint via an integrated fastfeedback internal backlight sensor, which continuously monitors and realigns these settings to maintain the factory calibration throughout the life of the monitor. This helps to reduce* the need for frequent manual or network-based calibration.

X-Light technology also alleviates the color shift of the backlight to the yellow spectrum, a common issue with competitive models. This benefit is achieved by keeping the whitepoint constant, while white light intensity remains the same. These two features are the basis for excellent diagnostic quality.


Monitor matching and reduced downtime. Another aspect of X-Light technology is its effective and efficient use in larger radiology departments and their IT infrastructure in hospitals. X-Light allows for easy multiple monitor matching — whether matching the whitepoints of two monitors or many monitors, the capability is virtually limitless. X-Light helps reduce service costs as it enables a more flexible use of display equipment. Since all X-Light-equipped monitors are set for the same whitepoint during production, they can be exchanged at any time within one location and swapped in existing dual- or multi-screen configurations. The commonly used method of manually screening displays for the appropriate matching whitepoint value and the need for an on-hand supply of "properly matched" exchange displays becomes an issue of the past. NEC MultiSync MD Series monitors have the ability to be matched to any other of their kind.

Standalone calibration and monitor matching. NEC MultiSync MD Series monitors have been designed to maintain their factory calibration. However,

if calibration or matching (Fig. 2) to a replacement is required, the monitors can be configured using the built-in, standalone algorithm programmed into the monitor's firmware. This feature is useful for monitors that are not used in a setup where calibration can easily be performed with a PC and external software or non-supported operating systems.

GammaComp™ MD calibration and conformance software. Quality control and calibration/matching can also be accomplished using GammaComp MD, ensuring consistent image quality on a single display system. The software maintains the monitor's conformance to the DICOM

Figure 2.



Display matching in the field

standard, while providing an easy-to-use QA environment for medical imaging. Optionally, GammaComp MD Administrator provides centralized control and management of multiple display systems.

3061 shades of gray. The NEC MultiSync MD Series provides support for improved diagnostic accuracy with simultaneous reproduction of up to 1024 from a palette of 3061 possible grayscales. The programmable 10-bit gamma correction ensures precise and smooth grayscale tuning and better representation of just noticeable differences (JNDs). Up to 1024 grayscales can be depicted at the same time for filmless diagnosis.



Figure 3. Attachable protective shield




 $\textbf{Figure 4.} \ \textbf{Sentry}^{\text{DX}} \ \textbf{attachable remote network sensor}$

Of course, grayscale images based on the internationally accepted DICOM standard (256 grayscales) can be viewed as well. Fast, internal switching of the related look-up tables (LUTs) supports the display's usage for softcopy reading of images from different medical modalities.

SA Superfine TFT displays. The innovative SA Superfine TFT glass with 2- or 3-Megapixel resolution delivers a brilliant picture with a wide viewing angle of up to 176° and minimal off-angle color shift. Therefore, it provides the best conditions for several doctors doing consultation in front of the display.

Anti-glare and low reflection. The SA Superfine TFT glass utilizes round resin beads instead of inconsistently shaped silica gel within the glass substrates for advanced glare reduction. Additionally, a low reflection overcoat is applied on top of the bead layer, which yields a low, uniform reflection over the surface of the glass. This process also provides for a truer representation of black and a higher contrast ratio.



A typical NEC MultiSync MD Series display system configuration

- 1 Single or dual head displays
- 2 Display controller board
- 3 GammaComp MD quality control software
- 4 Calibrator (optional)

Ultra-thin-frame design. The

dynamic cabinet design has become a trademark of NEC Display Solutions TFT displays. With a bezel width of only 16mm, they are perfectly suited to digital diagnostics systems in which two or more screens are positioned next to each other. The ultrathin frame minimizes eye distraction during side-by-side image viewing.

Mutiple video card support. NEC

MultiSync MD Series monitors work with 10-bit or 8-bit cards, however, 10-bit is recommended for optimum performance. 10-bit provides 1024 out of 3061 shades, while 8-bit provides only 256 out of 3061.

Internal power supply. The NEC MultiSync MD Series features an internal power supply. There is no

brick to kick, therefore, ground or desk clutter from a clunky power supply is eliminated.

Optional accessories. To add further performance and convenience, two accessories are available for the NEC MultiSync MD Series.


Protective shield. As these displays are often used by diagnostic professionals for illustration purposes, it is common for them to encounter frequent contact from fingers, pointers and pens. To help protect against LCD screen damage while maintaining the display's optimal image quality, an optional protective shield (Fig. 3) is available. This lightweight accessory, which is quickly and easily attached and removed, lengthens the life of the screen by fending off scratches caused by surface contact. The 100% transparent shield features anti-reflective coating on both sides, anti-fingerprint coating on the outside layer, a hard coat for ruggedness and excellent permanent marker resistance.

Sentry^{DX} remote network sensor. For non-assisted conformance, calibration and reporting functions, an optional, retractable remote sensor (Fig. 4) is available for NEC MultiSync MD Series displays. This component, which also can be used with non-MD Series monitors, is capable of measuring monitor brightness, whitepoint and ambient lighting, while performing a quality check for DICOM conformance and contrast response.

* Each end-user environment is unique, therefore NEC Display Solutions recommends that the user follows the quality control and calibration guidelines set forth by the AAPM TG-18.



When combined with an NEC Professional Series display (left), NEC MultiSync MD Series monitors can help provide an ideal PACS configuration for your facility.



Native Resolution Pixel Pitch 0.27mm 0.21mm Brightness at Native (typical) Contrast Ratio (typical) 900 cd/m2 (263 fl) max, 400 cd/m2 (117 fl) calibrated 700 cd/m2 (204 fl) max, 400 cd/m2 (117 fl) calibrated 700:1 700:1 Response Time (typical) 35ms Viewing Angle (typical) (up/down/left/right) 88°/88°/88°/88° 88°/88°/88°/88° Grayscale Tone 10-bit: 1024 shades of gray from a pallet of 3061 10-bit: 1024 shades of gray from a pallet of 3061 **Input Connectors** DVI-D and VGA Selectable Gamma DICOM, log-linear, 2.2, 1.8 and programmable DICOM, log-linear, 2.2, 1.8 and programmable White Point (color temperature) CB - P104 simulation, x=.292, y=.321, Approx. 8000 K CB - P104 simulation, x=.292, y=.321, Approx. 8000 K BB - P45 simulation, x=.265, y=.320, Approx. 14000 K BB - P45 simulation, x=.265, y=.320, Approx. 14000 K Sensor Calibration Brightness/White Point/External Yes/Yes/Yes Yes/Yes/Yes Power Supply (internal) Yes Yes Power Consumption (typical) (120V) 65W 65W **Power Savings** <3W <3W Tilt, Swivel Stand 25° up/5° down, 340° 25° up/5° down, 340° Height Adjustable Stand Pivot Enabled Stand 46mm 46mm Yes **VESA Mounting** 100 x 100mm 100 x 100mm **USB** Input For external colorimeter (calibration mode) For external colorimeter (calibration mode) Dimensions (WxHxD) 18.4 x 18.1 x 7.9 in. (Landscape) 18.4 x 18.1 x 7.9 in. (Landscape) 14.1 x 20.2 x 7.9 in. (Portrait) 14.1 x 20.2 x 7.9 in. (Portrait)

UL2601/EN60601-1/EC601, FCC part 15 class B, CE/MDD, PCT,

C-tick, PCBC/B Mark, PSB, EnergyStar, GEEA Energy Label, DIN

NOTE 1: American College of Radiology (ACR) recommends minimum of 70-ft. lamberts (fl)
NOTE 2: NEC MDview 19" and 21" color LCD monitors with UL60601-1 are available via special order.

16mm 25.8 lbs.

6868-57, FDA-510K

Real Vision, Tech Source Sentry^{DX} (MD-N2M5B)

1 year, parts and labor

EW1-MD3MP21/EW2-MD3MP21

MD21PS-BK

MDM10B-2MP

Display video card	MDM10B-2MP	MDM10B-3MP
Shipping Weight Shipping Dimensions (W x H x D) Graphics Chip Memory Type Memory Size Main RAMDAC Secondary RAMDAC Card Type Form Factor Maximum Resolution Connectors Certifications Limited Warranty	0.66 lbs. /0.3 kg 13.5 x 9.4 x 3 in. /344 x 239 x 76mm Matrox 512-bit GPU DDR (Double Data Rate) SDRAM 256 MB 400 MHz 400 MHz AGP 2x and 4x compatible or PCI 64-bit, 66 MHz ATX 1920 x 1200 (digital/analog) 2 x DVI-I Class A: FCC, CE, CSA, VCCI 1 year	0.66 lbs. /0.3 kg 13.5 x 9.4 x 3 in./ 344 x 239 x 76mm Matrox 512-bit GPU DDR (Double Data Rate) SDRAM 256 MB 400 MHz 400 MHz PCI 64-bit, 66 MHz ATX 2048 x 1536 (digital) 2 x DVI-I Class A: FCC, CE, CSA, VCCI 1 Year



Bezel Width

Net Weight

Regulatory Standards

Also Recommended

Limited Warranty

Protection Glass (optional)

Display Video Cards (preferred)

Remote Network Sensor (optional)

Extended Warranty (addtl. 1 yr./addtl. 2 yr.)

Model

LCD Viewable Image Size LCD Module Technology Active Screen Area



MultiSync is a registered trademark, and GammaComp, NaViSet and X-Light are trademarks of NEC Display Solutions. All other brand or product names are trademarks or registered trademarks of their respective holders. Product specifications subject to change.

@2005 NEC Display Solutions of America, Inc. All rights reserved.

16mm

25.8 lbs.

6868-57, FDA-510K

Real Vision, Tech Source

Sentry^{DX} (MD-N2M5B)

1 year, parts and labor

EW1-MD2MP21/EW2-MD2MP21

MD21PS-BK

MDM10B-2MP

NEC Display Solutions

500 Park Boulevard, Suite 1100 Itasca, IL 60143 866-NEC-MORE www.necdisplay.com



UL2601/EN60601-1/EC601, FCC part 15 class B, CE/MDD, PCT,

C-tick, PCBC/B Mark, PSB, EnergyStar, GEEA Energy Label, DIN